

PRIMASOFT[®]


PG4-1331

PROCEDURE EAR LOOP MASK
WHITE



PACK FACTORS
50/BOX; 500/CASE

1 LEVEL 1 BARRIER	
BFE (%)	≥ 99
PFE (%)	≥ 99
DELTA P Avg. (mmH2O/cm2)	≤2.7
FLUID RESISTANCE	80 mm Hg



MAIN USES:

- Protects both the patient and clinician from fluids and particles generated during medical procedures outside of the operating room, including wound care, line changes, dental work and more.
- Protects the wearer from contamination from droplets generated by coughing, talking, or sneezing of nearby persons.
- Also used for respiratory etiquette to contain droplets from symptomatic individuals.

FEATURES AND BENEFITS:

- PRIMASOFT™ masks are silky, smooth and irritant-free, ideal for those with sensitivities to other masks
- Full-width, aluminum alloy nose piece will not kink and ensures a secure fit over the nose and across the cheeks
- Soft, stretchy ear loops comfortably fit a wide variety of face sizes
- Suitable for respiratory etiquette and procedures where there is minimal risk of splash
- All seams are ultrasonically welded